


**AIM: To determine the calcium and magnesium hardness in a given water sample**

**Theory:** The property of water which restricts or checks the lather formation with soap is called hardness. In other words, the presence of multivalent cations, mostly calcium and magnesium ions, in water is referred to as hardness of water. Hardness is of two types: Temporary or carbonate hardness which can be removed by boiling and Permanent or Noncarbonate hardness which cannot be removed by boiling. The hardness is usually expressed in parts of  $\text{CaCO}_3$  equivalent or calcium and magnesium salts per million parts of water i.e in ppm.


The hardness of water can be determined by complexometric titration using Ethylene diammine tetra acetic acid (EDTA). EDTA in the form of its di-sodium salt forms complex with  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions of water sample. When Eriochrome Black T (EBT) indicator is added to the hard water at pH around 9–10, it gives wine red colored unstable complex with  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions of water sample. When this wine red colored complex is titrated against EDTA solution of known strength the  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions forms stable metal complex with EDTA and color changes from wine red to blue (color of EBT indicator) at the end point.




So titration at pH about 9–10 using EBT indicator gives the total amount of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions in the water sample.

Titration against EDTA at pH around 12.5 gives the hardness due to  $\text{Ca}^{2+}$  only. A pH of about 12.5 required for this titration can be obtained by adding diethyl amine base with 3-4 drops of calcon indicator or NaOH base with murexide indicator. At this high pH, the  $\text{Mg}^{2+}$  ion is quantitatively precipitated as  $\text{Mg}(\text{OH})_2$  and  $\text{Ca}^{2+}$  ion alone can be estimated by complexometric method using EDTA. At the end point color changes from pink to pure blue.


**Structure of di-sodium salt of EDTA:**



**Structure of EBT Indicator:**



**Structure of Calcon Indicator:**



**Apparatus:** Conical flask, burette, pipette, and measuring cylinder.

**Chemicals:** EDTA,  $\text{Zn}(\text{OAc})_2$ ,  $\text{NH}_4\text{Cl}$ ,  $\text{NH}_4\text{OH}$ , Eriochrome Black T (EBT), Calcon,  $\text{NaOH}$  and ethyl alcohol.

**Procedure:**

1. Primary standard  $\text{Zn}(\text{OAc})_2$  solution was provided.
2. Secondary standard EDTA solution and  $\text{NH}_4\text{Cl}$ – $\text{NH}_4\text{OH}$  buffer solution were provided.
3. **Standardization of EDTA solution:** Pipette out 10 mL of  $\text{Zn}(\text{OAc})_2$  solution in 100 mL conical flask, add 5 mL of  $\text{NH}_4\text{OH}$  solution (to neutralize the residual acid in the solution), 5 mL of buffer solution and 4–5 drops of EBT indicator. The solution becomes wine red color. Titrate with EDTA solution from the burette till the colour changes to clear blue. Repeat three times and find out the strength of EDTA from known strength of  $\text{Zn}(\text{OAc})_2$  solution.

4. **Estimation of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  in water:** Pipette out 20 mL of hard water sample in a 100 mL conical flask, add 5 mL of  $\text{NH}_4\text{OH}$  solution (to neutralize the residual acid in the solution), 5 mL of buffer solution and 4–5 drops of EBT indicator. Titrate with standard EDTA solution till wine red color changes to clear blue. Repeat three times.
5. **Estimation of  $\text{Ca}^{2+}$  in presence of  $\text{Mg}^{2+}$  in water:** Pipette out 20 mL of water sample in a 100 mL conical flask, add 3 mL of diethyl amine/ $\text{NaOH}$  solution. Shake the solution thoroughly to precipitate all  $\text{Mg}^{+2}$  ions as  $\text{Mg}(\text{OH})_2$ . Add 4-5drops of calcon indicator and titrate with standard EDTA solution till the pink colour of the solution changes to clear blue. Repeat three times.

### **Observations and Calculations:**

**Table1:- Standardization of EDTA solution**

Strength of  $\text{Zn}(\text{OAc})_2$  solution ( $S_1$ ) = 0.1(M)

| Entry | Volume of $\text{Zn}(\text{OAc})_2$ solution mL ( $V_1$ ) | Burette Reading (mL) |       | Volume of EDTA required (mL) | Mean volume of EDTA required ( $V_2$ ) (mL) |
|-------|---|----------------------|-------|------------------------------|---|
|       |   | Initial              | Final |                              |   |
| 1     | 10  |                      |       |                              |   |
| 2     | 10  |                      |       |                              |   |
| 3     | 10  |                      |       |                              |   |

$$V_1S_1 = V_2S_2$$

$$\therefore \text{Strength of EDTA solution } (S_2) = (V_1S_1 / V_2) \text{ (M)} = p(\text{M})$$

**Table 2:- Estimation of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  in given water sample**

| Entry | Volume of water sample taken (mL) | Burette Reading (mL) |       | Volume of EDTA required (mL) | Mean volume of EDTA required ( $V_m$ ) (mL) |
|-------|-----------------------------------|----------------------|-------|------------------------------|---|
|       |                                   | Initial              | Final |                              |   |
| 1     | 20                                |                      |       |                              |   |
| 2     | 20                                |                      |       |                              |   |
| 3     | 20                                |                      |       |                              |   |

EDTA forms 1:1 complex with  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$

1000 mL of 1(M) EDTA solution = 1 mole of  $\text{CaCO}_3$  = 100 g of  $\text{CaCO}_3$

1 mL of 1(M) EDTA solution = [100 x (1/1000)] g of  $\text{CaCO}_3$

$V_m$  mL of p(M) EDTA solution = [100 x (1/1000) x p x  $V_m$ ] g of  $\text{CaCO}_3$

= A g (say)

20 mL of water sample contain A g of  $\text{CaCO}_3$

1000 mL of water sample contain  $[(A/20) \times 1000]$  g of  $\text{CaCO}_3$

**Table 3:- Estimation of  $\text{Ca}^{2+}$  in given water sample**

| Entry | Volume of water sample taken (mL) | Burette Reading (mL) |       | Volume of EDTA required (mL) | Mean volume of EDTA required ( $V_n$ ) (mL) |
|-------|-----------------------------------|----------------------|-------|------------------------------|---|
|       |                                   | Initial              | Final |                              |   |
| 1     | 20                                |                      |       |                              |   |
| 2     | 20                                |                      |       |                              |   |
| 3     | 20                                |                      |       |                              |   |

1000 mL of 1(M) EDTA solution = 1 mole of  $\text{CaCO}_3$  = 100 g of  $\text{CaCO}_3$

1 mL of 1(M) EDTA solution =  $[100 \times (1/1000)]$  g of  $\text{CaCO}_3$

$(V_m - V_n)$  mL of p(M) EDTA solution =  $[100 \times (1/1000) \times p \times (V_m - V_n)]$  g of  $\text{CaCO}_3$   
= B g (say)

20 mL of water sample contain B g of  $\text{CaCO}_3$

1000 mL of water sample contain  $[(B/20) \times 1000]$  g of  $\text{CaCO}_3$

**Conclusion:** Total hardness of water sample w.r.t  $\text{CaCO}_3$  =  $[(A/20) \times 1000]$  g/L

Magnesium hardness of water w.r.t  $\text{CaCO}_3$  =  $[(B/20) \times 1000]$  g/L

Calcium hardness of water w.r.t  $\text{CaCO}_3$  =  $[(A/20) \times 1000] - [(B/20) \times 1000]$  g/L

### **Discussion:**

To determine the temporary and permanent hardness separately in the water sample, same EDTA method can be used first to determine the total hardness ( $H_1$ ). Then boil the sample water gently for about one hour, cool it and filter it. Take this sample of water and estimate the permanent hardness ( $H_2$ ) by the same EDTA method. As it is known that, temporary hardness can be removed just by boiling. The difference ( $H_1 - H_2$ ) is the temporary hardness in the given water sample.

### **Precautions:**

- All the glass apparatus should be washed thoroughly with distilled water before use.
- The burette and pipette should be rinsed with solution to be taken in it.
- There should not be any leakage in the burette.
- All the solution should be freshly prepared.
- Same amount of indicator should be added each time.

- vi) pH of the solution should be maintained during titration.
- vii) Shaking should be proper during titration.
- viii) The titration flask should be placed on white paper or board to identify properly the colour change at the end point.